

产品名称: Annexin V Rabbit Monoclonal Antibody
产品货号: AMRe85204



产品概述 (Summary)

产品名称 (Production Name)	Annexin V Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Annexin V
别名 (Alternative Names)	PP4; ANX5; ENX2; RPRGL3.
基因 ID (Gene ID)	308.0
蛋白 ID (SwissProt ID)	P08758.A synthetic peptide of human Annexin V

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 36 kDa

研究背景 (Background)

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is


产品名称: Annexin V Rabbit Monoclonal Antibody
产品货号: AMRe85204




involved in the blood coagulation cascade.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Annexin V in K562, C6, 3T3, HeLa lysates using Annexin V antibody.



Immunocytochemistry analysis of Annexin V(green) in HeLa using Annexin V antibody, and DAPI(blue)

注意事项 (Note)

For research use only .