

产品名称: USP11 Rabbit Monoclonal Antibody
产品货号: AMRe85213



产品概述 (Summary)

产品名称 (Production Name)	USP11 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP11
别名 (Alternative Names)	UHX1; USP11
基因 ID (Gene ID)	8237.0
蛋白 ID (SwissProt ID)	P51784.A synthetic peptide of human USP11

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 110 kDa; Observed MW: 110 kDa

研究背景 (Background)

Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Inhibits the degradation of


产品名称: USP11 Rabbit Monoclonal Antibody
产品货号: AMRe85213



target proteins by the proteasome. Plays a role in the regulation of pathways leading to NF-kappa-B activation. Plays a role in the regulation of DNA repair after double-stranded DNA breaks.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of USP11 in Jurkat, rat Brain, HeLa lysates using USP11 antibody.

注意事项 (Note)

For research use only .