

产品名称: NEDD8 Rabbit Monoclonal Antibody
产品货号: AMRe85218



产品概述 (Summary)

产品名称 (Production Name)	NEDD8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEDD8
别名 (Alternative Names)	NEDD8; NEDD8; Neddylin; Neural precursor cell expressed developmentally down-regulated protein 8; NEDD-8; Ubiquitin-like protein Nedd8
基因 ID (Gene ID)	4738.0
蛋白 ID (SwissProt ID)	Q15843.A synthetic peptide of human NEDD8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 9 kDa; Observed MW: 9 kDa

研究背景 (Background)


产品名称: NEDD8 Rabbit Monoclonal Antibody
产品货号: AMRe85218



Neural precursor cell-expressed developmentally downregulated protein 8 (NEDD8), also known as Rub1 (related to ubiquitin 1) in plants and yeast, is a member of the ubiquitin-like protein family. The covalent attachment of NEDD8 to target proteins, termed neddylation, is a reversible, multi-step process analogous to ubiquitination.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of NEDD8 in Jurkat lysates using NEDD8 antibody.



Immunocytochemistry analysis of NEDD8 (green) in HeLa using NEDD8 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NEDD8 antibody. High-pressure and

产品名称: NEDD8 Rabbit Monoclonal Antibody
产品货号: AMRe85218



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .