

产品名称: CD8 alpha Rabbit Monoclonal Antibody
产品货号: AMRe85220



产品概述 (Summary)

产品名称 (Production Name)	CD8 alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD8 alpha
别名 (Alternative Names)	CD8A; MAL; T-cell surface glycoprotein CD8 alpha chain; T-lymphocyte differentiation antigen T8/Leu-2; CD antigen CD8a
基因 ID (Gene ID)	925.0
蛋白 ID (SwissProt ID)	P01732.A synthetic peptide of human CD8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 29 kDa

研究背景 (Background)


产品名称: CD8 alpha Rabbit Monoclonal Antibody
产品货号: AMRe85220




Identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing. CD8 alpha chains binds to class I MHC molecules alpha-3 domains.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of CD8 in THP-1 lysates using CD8 alpha antibody.



Western blot analysis of CD8 in K562, CHO-K1, HeLa lysates using CD8 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: CD8 alpha Rabbit Monoclonal Antibody
产品货号: AMRe85220



注意事项 (Note)

For research use only .