

产品名称: BNIP3 Rabbit Monoclonal Antibody  
产品货号: AMRe85223



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | BNIP3 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB, ICC                                |
| 种属反应性 (Reactivity)     | Human, Rat                             |

## 产品性能 (Performance)

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                   |
| 形式 (Form)           | Liquid                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.     |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                           |
|------------------------|-------------------------------------------|
| 基因名 (Gene Name)        | BNIP3                                     |
| 别名 (Alternative Names) | BNIP3; NIP3                               |
| 基因 ID (Gene ID)        | 664.0                                     |
| 蛋白 ID (SwissProt ID)   | Q12983.A synthetic peptide of human BNIP3 |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000, ICC 1:50-1:200            |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 22 kDa; Observed MW: 28 kDa |

## 研究背景 (Background)

Apoptosis-inducing protein that, which can overcome BCL2 suppression. May play a role in repartitioning calcium between

产品名称: BNIP3 Rabbit Monoclonal Antibody  
产品货号: AMRe85223




the two major intracellular calcium stores in association with BCL2.


## 研究领域 (Research Area)

Apoptosis, mTOR signaling pathway


## 图片 (Image Data)



Western blot analysis of BNIP3 in SH-SY5Y lysates using BNIP3 antibody.



Western blot analysis of BNIP3 in rat Brain lysates using BNIP3 antibody.



Immunocytochemistry analysis of BNIP3(green) in HeLa using BNIP3 antibody, and DAPI(blue)

**产品名称: BNIP3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe85223**



---

**注意事项 (Note)**

For research use only .