

产品名称: NLRP3 Rabbit Monoclonal Antibody
产品货号: AMRe85229



产品概述 (Summary)

产品名称 (Production Name)	NLRP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NLRP3
别名 (Alternative Names)	FCU; MWS; FCAS; Cias1; Mmig1; NLRP3; Pypaf1; Ail/AVP; AGTAVPRL; Cryopyrin
基因 ID (Gene ID)	114548.0
蛋白 ID (SwissProt ID)	Q96P20.Recombinant protein of human NLRP3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 118 kDa; Observed MW: 118 kDa

研究背景 (Background)


产品名称: NLRP3 Rabbit Monoclonal Antibody
产品货号: AMRe85229



As the sensor component of the NLRP3 inflammasome, plays a crucial role in innate immunity and inflammation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NLRP3 in THP-1 lysates using NLRP3 antibody.

注意事项 (Note)

For research use only .