

产品名称: c-Jun Rabbit Monoclonal Antibody
产品货号: AMRe85237



产品概述 (Summary)

产品名称 (Production Name)	c-Jun Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	c-Jun
别名 (Alternative Names)	Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39
基因 ID (Gene ID)	3725.0
蛋白 ID (SwissProt ID)	P05412.A synthetic peptide of human c-Jun

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 48,43 kDa

研究背景 (Background)

产品名称: c-Jun Rabbit Monoclonal Antibody
产品货号: AMRe85237




c-jun is a transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA[CG]TCA-3'. Promotes activity of NR5A1 when phosphorylated by HIPK3 leading to increased steroidogenic gene expression upon cAMP signaling pathway stimulation.


研究领域 (Research Area)

TGF-beta signaling pathway, MAPK signaling pathway


图片 (Image Data)



Western blot analysis of c-Jun in 293, 293T lysates using c-Jun antibody.



Western blot analysis of c-Jun in 3T3 lysates using c-Jun antibody



Immunocytochemistry analysis of c-Jun (green) in 293T using c-Jun antibody, and DAPI (blue).

产品名称: c-Jun Rabbit Monoclonal Antibody
产品货号: AMRe85237



注意事项 (Note)

For research use only .