

产品名称: 14-3-3 beta/alpha Rabbit Monoclonal Antibody
产品货号: AMRe85242



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 beta/alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	14-3-3 protein beta/alpha
别名 (Alternative Names)	YWHAB; 14-3-3 protein beta/alpha; Protein 1054; Protein kinase C inhibitor protein 1; KCIP-1
基因 ID (Gene ID)	7529.0
蛋白 ID (SwissProt ID)	P31946.A synthetic peptide of human 14-3-3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

研究背景 (Background)


产品名称: 14-3-3 beta/alpha Rabbit Monoclonal Antibody
产品货号: AMRe85242




Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner. Negative regulator of osteogenesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 41701 in K562, rat Brain, C6, 3T3, HeLa lysates using 14-3-3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using 14-3-3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .