

产品名称: Activin Receptor Type IA Rabbit Monoclonal Antibody
产品货号: AMRe85254



产品概述 (Summary)

产品名称 (Production Name)	Activin Receptor Type IA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Activin Receptor Type IA
别名 (Alternative Names)	ACTRI; Acvr1; ACVR1A; ACVRLK2; ALK2; FOP; SKR1; TSRI
基因 ID (Gene ID)	90.0
蛋白 ID (SwissProt ID)	Q04771.A synthetic peptide of human Activin Receptor Type IA

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 57 kDa

研究背景 (Background)


产品名称: Activin Receptor Type IA Rabbit Monoclonal Antibody
产品货号: AMRe85254




On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Activin Receptor Type IA in mouse heart, mouse brain lysates using Activin Receptor Type IA antibody.



Western blot analysis of Activin Receptor Type IA in HeLa lysates using Activin Receptor Type IA antibody.

注意事项 (Note)

For research use only .