

产品名称: ALKBH1 Rabbit Monoclonal Antibody
产品货号: AMRe85272



产品概述 (Summary)

产品名称 (Production Name)	ALKBH1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALKBH1
别名 (Alternative Names)	ABH; ABH1; alkB; ALKBH; ALKBH1; hABH
基因 ID (Gene ID)	8846.0
蛋白 ID (SwissProt ID)	Q13686.Recombinant protein of human ALKBH1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 44 kDa

研究背景 (Background)

Dioxygenase that repairs alkylated single-stranded DNA and RNA containing 3-methylcytosine by oxidative demethylation.


产品名称: ALKBH1 Rabbit Monoclonal Antibody
产品货号: AMRe85272




Requires molecular oxygen, alpha-ketoglutarate and iron.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ALKBH1 in 3T3 lysates using ALKBH1 antibody.



Western blot analysis of ALKBH1 in Jurkat, C6, HeLa lysates using ALKBH1 antibody.

注意事项 (Note)

For research use only .