

产品名称: alpha Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe85279



产品概述 (Summary)

产品名称 (Production Name)	alpha Tubulin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	alpha Tubulin
别名 (Alternative Names)	ALS22; TUBA1; H2-ALPHA
基因 ID (Gene ID)	7277.0
蛋白 ID (SwissProt ID)	P68366.A synthetic peptide of human alpha Tubulin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

研究背景 (Background)

TUBA1A Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the


产品名称: alpha Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe85279




beta chain and one at a non-exchangeable site on the alpha-chain.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of alpha Tubulin in mouse brain, 3T3 lysates using alpha Tubulin antibody.




Western blot analysis of alpha Tubulin in HeLa, CHO-K1, rat Brain lysates using alpha Tubulin antibody.



Immunocytochemistry analysis of alpha Tubulin (green) in HeLa cells using alpha Tubulin antibody and DAPI (blue).

产品名称: alpha Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe85279



Immunohistochemistry analysis of paraffin-embedded Human tonsil using alpha Tubulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .