

产品名称: Aly Rabbit Monoclonal Antibody
产品货号: AMRe85283



产品概述 (Summary)

产品名称 (Production Name)	Aly Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Aly ALYREF; ALY; BEF; THOC4; THO complex subunit 4; Tho4; Ally of AML-1 and LEF-1; Aly/REF export factor; Transcriptional coactivator Aly/REF; bZIP-enhancing factor BEF
别名 (Alternative Names)	
基因 ID (Gene ID)	10189.0
蛋白 ID (SwissProt ID)	Q86V81.A synthetic peptide of human Aly

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 27 kDa; Observed MW: 27 kDa

产品名称: Aly Rabbit Monoclonal Antibody
产品货号: AMRe85283




研究背景 (Background)


Component of the THO subcomplex of the TREX complex. The TREX complex specifically associates with spliced mRNA and not with unspliced pre-mRNA. It is recruited to spliced mRNAs by a transcription-independent mechanism. Binds to mRNA upstream of the exon-junction complex (EJC) and is recruited in a splicing- and cap-dependent manner to a region near the 5' end of the mRNA where it functions in mRNA export.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Aly/Ref in rat brain, rat heart lysates using Aly antibody.



Western blot analysis of Aly/Ref in HepG2, mouse liver lysates using Aly/Ref antibody.

注意事项 (Note)

For research use only .