

产品名称: Annexin A2 Rabbit Monoclonal Antibody
产品货号: AMRe85288



产品概述 (Summary)

产品名称 (Production Name)	Annexin A2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Annexin A2 ANXA2; ANX2; ANX2L4; CAL1H; LPC2D; Annexin A2; Annexin II; Annexin-2;
别名 (Alternative Names)	Calpactin I heavy chain; Calpactin-1 heavy chain; Chromobindin-8; Lipocortin II; Placental anticoagulant protein IV; PAP-IV; Protein I; p36
基因 ID (Gene ID)	302.0
蛋白 ID (SwissProt ID)	P07355.A synthetic peptide of human Annexin A2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 39 kDa

产品名称: Annexin A2 Rabbit Monoclonal Antibody
产品货号: AMRe85288




研究背景 (Background)


Originally identified as a protein inhibitor of phospholipase A2, annexin A2 has subsequently been shown to interact with an array of protein and non-protein partners, including F-actin, spectrin, SNARE complexes, RNA, and virus particles.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Annexin A2 in C6, 3T3, HeLa lysates using Annexin A2 antibody.



Immunocytochemistry analysis of Annexin A2(green) in HeLa using Annexin A2 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .