

产品名称: Transcription Factor AP 2 gamma Rabbit
Monoclonal Antibody
产品货号: AMRe85292



产品概述 (Summary)

产品名称 (Production Name)	Transcription Factor AP 2 gamma Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Transcription Factor AP 2 gamma
别名 (Alternative Names)	ERF1; TFAP2G; hAP-2g; AP2-GAMMA
基因 ID (Gene ID)	7022.0
蛋白 ID (SwissProt ID)	Q92754. Recombinant protein of human AP2 gamma

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa; Observed MW: 49 kDa


研究背景 (Background)

产品名称: Transcription Factor AP 2 gamma Rabbit Monoclonal Antibody
产品货号: AMRe85292


The protein encoded by this gene is a sequence-specific DNA-binding transcription factor involved in the activation of several developmental genes. The encoded protein can act as either a homodimer or heterodimer with other family members and is induced during retinoic acid-mediated differentiation. It plays a role in the development of the eyes, face, body wall, limbs, and neural tube.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TFAP2C/AP2 gamma in Jurkat, HeLa lysates using Transcription Factor AP 2 gamma antibody.



Immunocytochemistry analysis of TFAP2C/AP2 gamma (green) in MCF-7 using TFAP2C/AP2 gamma antibody, and DAPI (blue)

注意事项 (Note)

For research use only .