

产品名称: AP3M1 Rabbit Monoclonal Antibody
产品货号: AMRe85293



产品概述 (Summary)

产品名称 (Production Name)	AP3M1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AP3M1
别名 (Alternative Names)	AP-3 adaptor complex mu3A subunit; Adaptor-related protein complex 3 subunit mu-1; Mu-adaptin 3A; Mu3A-adaptin
基因 ID (Gene ID)	26985.0
蛋白 ID (SwissProt ID)	Q9Y2T2.A synthetic peptide of human AP3M1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)


产品名称: AP3M1 Rabbit Monoclonal Antibody
产品货号: AMRe85293



Part of the AP-3 complex, an adaptor-related complex which is not clathrin-associated. The complex is associated with the Golgi region as well as more peripheral structures. It facilitates the budding of vesicles from the Golgi membrane and may be directly involved in trafficking to lysosomes. In concert with the BLOC-1 complex, AP-3 is required to target cargos into vesicles assembled at cell bodies for delivery into neurites and nerve terminals.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of AP3M1 in K562, C6, 3T3, HeLa lysates using AP3M1 antibody.

注意事项 (Note)

For research use only .