

产品名称: ATG12 Rabbit Monoclonal Antibody
产品货号: AMRe85295



产品概述 (Summary)

产品名称 (Production Name)	ATG12 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATG12
别名 (Alternative Names)	APG12-like; APG12L; ATG12; Autophagy 12; Autophagy-related protein 12; FBR93; HAPG12; Ubiquitin-like protein ATG12
基因 ID (Gene ID)	9140.0
蛋白 ID (SwissProt ID)	O94817.A synthetic peptide of human Apg12

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 55 kDa

研究背景 (Background)

产品名称: ATG12 Rabbit Monoclonal Antibody
产品货号: AMRe85295




Ubiquitin-like protein required for autophagy. Conjugated to ATG3 and ATG5.

研究领域 (Research Area)

Autophagy

图片 (Image Data)



Western blot analysis of Apg12 in K562, 3T3, HeLa lysates using ATG12 antibody.

注意事项 (Note)

For research use only .