

产品名称: Actin Related Protein 3 Rabbit Monoclonal Antibody
产品货号: AMRe85307



产品概述 (Summary)

产品名称 (Production Name)	Actin Related Protein 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Actin Related Protein 3
别名 (Alternative Names)	ACTR3; ARP3; Actin-related protein 3; Actin-like protein 3
基因 ID (Gene ID)	10096.0
蛋白 ID (SwissProt ID)	P61158.A synthetic peptide of human Arp3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)


产品名称: Actin Related Protein 3 Rabbit Monoclonal Antibody
产品货号: AMRe85307




The Arp2/3 complex involved in Golgi polarization in NIH3T3 cells. In a different model it has been shown that the Arp2/3 complex is necessary for neutrophil chemotaxis and phagocytosis.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Arp3 in K562, rat Brain, C6, 3T3, HeLa lysates using Arp3 antibody.



Immunocytochemistry analysis of Arp3(green) in HeLa cells using Arp3 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded mouse colon using Actin Related Protein 3 antibody. High-pressure

**产品名称: Actin Related Protein 3 Rabbit Monoclonal
Antibody**
产品货号: AMRe85307



and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .