

产品名称: Aryl Hydrocarbon Receptor Rabbit Monoclonal Antibody
产品货号: AMRe85311



产品概述 (Summary)

产品名称 (Production Name)	Aryl Hydrocarbon Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Aryl Hydrocarbon Receptor
别名 (Alternative Names)	Ah receptor; AhR; Class E basic helix-loop-helix protein 76; bHLHe76; AHR
基因 ID (Gene ID)	196.0
蛋白 ID (SwissProt ID)	P35869.A synthetic peptide of human Aryl hydrocarbon Receptor

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa; Observed MW: 100 kDa


研究背景 (Background)

产品名称: Aryl Hydrocarbon Receptor Rabbit Monoclonal Antibody
产品货号: AMRe85311


The aryl hydrocarbon receptor (AhR) is a ligand activated transcription factor involved in xenobiotic metabolism, cell cycle regulation, and development in response to both endogenous and environmental signals. Involved in cell-cycle regulation. Likely to play an important role in the development and maturation of many tissues.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Aryl hydrocarbon Receptor in HeLa lysates using Aryl Hydrocarbon Receptor antibody.



Immunocytochemistry analysis of Aryl hydrocarbon Receptor (green) in HeLa using Aryl hydrocarbon Receptor antibody, and DAPI (blue)

注意事项 (Note)

For research use only .