

产品名称: ASK1 Rabbit Monoclonal Antibody
产品货号: AMRe85314



产品概述 (Summary)

产品名称 (Production Name)	ASK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASK1 MAP3K5; ASK1; MAPKKK5; MEKK5; Mitogen-activated protein kinase kinase
别名 (Alternative Names)	kinase 5; Apoptosis signal-regulating kinase 1; ASK-1; MAPK/ERK kinase kinase 5; MEK kinase 5; MEKK 5
基因 ID (Gene ID)	4217.0
蛋白 ID (SwissProt ID)	Q99683.A synthetic peptide of human ASK1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 155 kDa; Observed MW: 155 kDa

产品名称: ASK1 Rabbit Monoclonal Antibody
产品货号: AMRe85314




研究背景 (Background)

Mitogen-activated protein kinase (MAPK) signaling cascades include MAPK or extracellular signal-regulated kinase (ERK), MAPK kinase (MKK or MEK), and MAPK kinase kinase (MAPKKK or MEKK). MAPKK kinase/MEKK phosphorylates and activates its downstream protein kinase, MAPK kinase/MEK, which in turn activates MAPK.

研究领域 (Research Area)

Apoptosis, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of ASK1 in HeLa lysates using ASK1 antibody.

注意事项 (Note)

For research use only .