

产品名称: Aspartate Aminotransferase Rabbit Monoclonal Antibody
产品货号: AMRe85315

产品概述 (Summary)

产品名称 (Production Name)	Aspartate Aminotransferase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Aspartate Aminotransferase
别名 (Alternative Names)	Aspartate aminotransferase; cytoplasmic; Glutamate oxaloacetate transaminase 1; Transaminase A
基因 ID (Gene ID)	2805.0
蛋白 ID (SwissProt ID)	P17174.A synthetic peptide of human Aspartate Aminotransferase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa; Observed MW: 46 kDa

产品名称: Aspartate Aminotransferase Rabbit Monoclonal Antibody
产品货号: AMRe85315




研究背景 (Background)


Plays a key role in amino acid metabolism.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Aspartate Aminotransferase in K562, rat Brain, 3T3, HeLa lysates using Aspartate Aminotransferase antibody.



Immunocytochemistry analysis of Aspartate Aminotransferase (green) in U87-MG using Aspartate Aminotransferase antibody, and DAPI (blue).

产品名称: Aspartate Aminotransferase Rabbit Monoclonal Antibody
产品货号: AMRe85315



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using Aspartate Aminotransferase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .