

产品名称: Axin 2 Rabbit Monoclonal Antibody
产品货号: AMRe85331



产品概述 (Summary)

产品名称 (Production Name)	Axin 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Axin 2
别名 (Alternative Names)	Axil; Axin like protein; Axin2; Axis inhibition protein 2; Conductin
基因 ID (Gene ID)	8313.0
蛋白 ID (SwissProt ID)	Q9Y2T1.A synthetic peptide of human Axin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa; Observed MW: 94 kDa

研究背景 (Background)

Axin1 (Axis inhibition protein 1) and Axin2 are multidomain scaffold proteins that negatively regulate Wnt signaling. Axin1

产品名称: Axin 2 Rabbit Monoclonal Antibody
产品货号: AMRe85331




interacts with APC, GSK-3 β , Dvl, and β -catenin and promotes the GSK-3 β -mediated phosphorylation and subsequent degradation of β -catenin .Down-regulates beta-catenin.

研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Western blot analysis of Axin 2 in K562, C6, 3T3, HeLa lysates using Axin 2 antibody.

注意事项 (Note)

For research use only .