

产品名称: AXL Rabbit Monoclonal Antibody
产品货号: AMRe85332



产品概述 (Summary)

产品名称 (Production Name)	AXL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AXL
别名 (Alternative Names)	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
基因 ID (Gene ID)	558.0
蛋白 ID (SwissProt ID)	P30530.Recombinant protein of human Axl

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 98 kDa; Observed MW: 138 kDa

研究背景 (Background)

The protein encoded by this gene is a member of the Tyro3-Axl-Mer (TAM) receptor tyrosine kinase subfamily. The


产品名称: AXL Rabbit Monoclonal Antibody
产品货号: AMRe85332



encoded protein possesses an extracellular domain which is composed of two immunoglobulin-like motifs at the N-terminal, followed by two fibronectin type-III motifs. It transduces signals from the extracellular matrix into the cytoplasm by binding to the vitamin K-dependent protein growth arrest-specific 6 (Gas6). This gene may be involved in several cellular functions including growth, migration, aggregation and anti-inflammation in multiple cell types. Alternative splicing results in multiple transcript variants of this gene.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of AXL Rabbit mAb in Hela lysates using AXL antibody.

注意事项 (Note)

For research use only .