

产品名称: Bad Rabbit Monoclonal Antibody
产品货号: AMRe85336



产品概述 (Summary)

产品名称 (Production Name)	Bad Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bad BAD;BBC6; BCL2L8; Bcl2 antagonist of cell death; BAD; Bcl-2-binding component 6; Bcl-2-like protein 8; Bcl2-L-8; Bcl-XL/Bcl-2-associated death promoter
别名 (Alternative Names)	
基因 ID (Gene ID)	572.0
蛋白 ID (SwissProt ID)	Q92934.A synthetic peptide of human Bad

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa; Observed MW: 23 kDa

产品名称: Bad Rabbit Monoclonal Antibody
产品货号: AMRe85336




研究背景 (Background)

Promotes cell death. Successfully competes for the binding to Bcl-X(L), Bcl-2 and Bcl-W, thereby affecting the level of heterodimerization of these proteins with BAX. Can reverse the death repressor activity of Bcl-X(L), but not that of Bcl-2 (By similarity). Appears to act as a link between growth factor receptor signaling and the apoptotic pathways.


研究领域 (Research Area)

Apoptosis, PI3K-Akt signaling pathway


图片 (Image Data)



Western blot analysis of Bad in rat Brain, C6, 3T3, HeLa lysates using Bad antibody.



Immunocytochemistry analysis of Bad(green) in HeLa using Bad antibody, and DAPI(blue)



产品名称: Bad Rabbit Monoclonal Antibody
产品货号: AMRe85336



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Bad antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .