

产品名称: BAG1 Rabbit Monoclonal Antibody
产品货号: AMRe85337



产品概述 (Summary)

产品名称 (Production Name)	BAG1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAG1
别名 (Alternative Names)	BAG1; HAP; BAG family molecular chaperone regulator 1; BAG-1; Bcl-2-associated athanogene 1
基因 ID (Gene ID)	573.0
蛋白 ID (SwissProt ID)	Q99933.A synthetic peptide of human Bag1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 52, 46, 33 kDa

研究背景 (Background)

产品名称: BAG1 Rabbit Monoclonal Antibody
产品货号: AMRe85337




Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Inhibits the pro-apoptotic function of PPP1R15A, and has anti-apoptotic activity. Markedly increases the anti-cell death function of BCL2 induced by various stimuli.


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of Bag1 in HeLa, MCF-7, U251MG, K562 lysates using BAG1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Bag1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .