

产品名称: BAT3 Rabbit Monoclonal Antibody
产品货号: AMRe85342



产品概述 (Summary)

产品名称 (Production Name)	BAT3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAT3
别名 (Alternative Names)	BAG6; BAT 3; Scythe
基因 ID (Gene ID)	7917.0
蛋白 ID (SwissProt ID)	P46379.A synthetic peptide of human BAT3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 119 kDa; Observed MW: 150 kDa

研究背景 (Background)

Chaperone that plays a key role in various processes such as apoptosis, insertion of tail-anchored (TA) membrane proteins

产品名称: BAT3 Rabbit Monoclonal Antibody
产品货号: AMRe85342




to the endoplasmic reticulum membrane and regulation of chromatin. Acts in part by regulating stability of proteins and their degradation by the proteasome. Participates in endoplasmic reticulum stress-induced apoptosis via its interaction with AIFM1/AIF by regulating AIFM1/AIF stability and preventing its degradation.


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of BAT3/BAG6 in K562, rat Brain, C6, 3T3 lysates using BAT3 antibody.



Immunocytochemistry analysis of BAT3/BAG6 (green) in MCF-7 using BAT3/BAG6 antibody, and DAPI (blue)

注意事项 (Note)

For research use only .