

**产品名称: Bcl10 Rabbit Monoclonal Antibody**  
**产品货号: AMRe85344**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Bcl10 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB                                     |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                               |
| 修饰 (Modification)   | Unmodified                                                                                 |
| 同种型 (Isotype)       | IgG                                                                                        |
| 克隆 (Clonality)      | Monoclonal                                                                                 |
| 形式 (Form)           | Liquid                                                                                     |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.   |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                      |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                         |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | Bcl10<br>BCL10; CIPER; CLAP; B-cell lymphoma/leukemia 10; B-cell CLL/lymphoma 10;                                                                                       |
| 别名 (Alternative Names) | Bcl-10; CARD-containing molecule enhancing NF-kappa-B; CARD-like apoptotic protein; hCLAP; CED-3/ICH-1 prodomain homologous E10-like regulator; CIPER; Cellular homolog |
| 基因 ID (Gene ID)        | 8915.0                                                                                                                                                                  |
| 蛋白 ID (SwissProt ID)   | O95999.A synthetic peptide of human Bcl10                                                                                                                               |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000                            |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 26 kDa; Observed MW: 26 kDa |

产品名称: Bcl10 Rabbit Monoclonal Antibody  
产品货号: AMRe85344




## 研究背景 (Background)

Promotes apoptosis, pro-caspase-9 maturation and activation of NF-kappa-B via NIK and IKK. May be an adapter protein between upstream TNFR1-TRADD-RIP complex and the downstream NIK-IKK-IKAP complex. Is a substrate for MALT1.

## 研究领域 (Research Area)

Apoptosis

## 图片 (Image Data)



Western blot analysis of Bcl10 in K562, C6, 3T3, HeLa lysates using Bcl10 antibody.

## 注意事项 (Note)

For research use only .