

产品名称: Bcl-XL Rabbit Monoclonal Antibody
产品货号: AMRe85345



产品概述 (Summary)

产品名称 (Production Name)	Bcl-XL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bcl-XL
别名 (Alternative Names)	BCL2L1; BCL2L; BCLX; Bcl-2-like protein 1; Bcl2-L-1; Apoptosis regulator Bcl-X
基因 ID (Gene ID)	598.0
蛋白 ID (SwissProt ID)	Q07817.Recombinant protein of human Bcl-XL

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 30 kDa

研究背景 (Background)

The protein encoded by this gene belongs to the BCL-2 protein family. BCL-2 family members form hetero- or homodimers

产品名称: Bcl-XL Rabbit Monoclonal Antibody
产品货号: AMRe85345




and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities.


研究领域 (Research Area)

Apoptosis, PI3K-Akt signaling pathway


图片 (Image Data)



Western blot analysis of Bcl-XL in rat Brain, C6, 3T3, HeLa lysates using Bcl-XL antibody.




Western blot analysis of Bcl-XL in K562, rat Brain, 3T3 lysates using Bcl-XL antibody



Immunocytochemistry analysis of Bcl-XL (green) in HT-1080 using Bcl-XL antibody, and DAPI (blue).

产品名称: Bcl-XL Rabbit Monoclonal Antibody
产品货号: AMRe85345



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Bcl-XL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .