

产品名称: BMP2 Rabbit Monoclonal Antibody
产品货号: AMRe85354



产品概述 (Summary)

产品名称 (Production Name)	BMP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BMP2
别名 (Alternative Names)	BMP2; BMP2A; Bone morphogenetic protein 2; BMP-2; Bone morphogenetic protein 2A; BMP-2A
基因 ID (Gene ID)	650.0
蛋白 ID (SwissProt ID)	P12643.Recombinant protein of human BMP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 30 kDa

研究背景 (Background)

产品名称: BMP2 Rabbit Monoclonal Antibody
产品货号: AMRe85354




Bone morphogenic proteins (BMPs) are members of the TGF β superfamily. BMPs are involved in the induction of cartilage and bone formation. Induces cartilage and bone formation.

研究领域 (Research Area)

TGF-beta signaling pathway

图片 (Image Data)



Western blot analysis of BMP2 in MCF-7, HeLa lysates using BMP2 antibody.

注意事项 (Note)

For research use only .