

产品名称: Carbonic Anhydrase 1 Rabbit Monoclonal Antibody
产品货号: AMRe85376



产品概述 (Summary)

产品名称 (Production Name)	Carbonic Anhydrase 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Carbonic Anhydrase 1
别名 (Alternative Names)	CA1; Carbonic anhydrase 1; Carbonate dehydratase I; Carbonic anhydrase B; CAB; Carbonic anhydrase I; CA-I
基因 ID (Gene ID)	759.0
蛋白 ID (SwissProt ID)	P00915.A synthetic peptide of human Carbonic Anhydrase I

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 29 kDa

产品名称: Carbonic Anhydrase 1 Rabbit Monoclonal Antibody
产品货号: AMRe85376




研究背景 (Background)


Reversible hydration of carbon dioxide. Can hydrates cyanamide to urea.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Carbonic Anhydrase 1/CA1 in TF1, mouse gastric lysates using Carbonic Anhydrase 1 antibody.



Western blot analysis of Carbonic Anhydrase 1/CA1 in rat Brain lysates using Carbonic Anhydrase 1/CA1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Carbonic Anhydrase I antibody. High-

产品名称: Carbonic Anhydrase 1 Rabbit Monoclonal Antibody
产品货号: AMRe85376



pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .