

产品名称: Integrin alpha L Rabbit Monoclonal Antibody
产品货号: AMRe85392



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha L Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin alpha L ITGAL; CD11A; Integrin alpha-L; CD11 antigen-like family member A;
别名 (Alternative Names)	Leukocyte adhesion glycoprotein LFA-1 alpha chain; LFA-1A; Leukocyte function-associated molecule 1 alpha chain; CD11a
基因 ID (Gene ID)	3683.0
蛋白 ID (SwissProt ID)	P20701.A synthetic peptide of human CD11a

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 129 kDa; Observed MW: 150 kDa

产品名称: Integrin alpha L Rabbit Monoclonal Antibody
产品货号: AMRe85392




研究背景 (Background)

Integrin alpha-L/beta-2 is a receptor for ICAM1, ICAM2, ICAM3 and ICAM4. It is involved in a variety of immune phenomena including leukocyte-endothelial cell interaction, cytotoxic T-cell mediated killing, and antibody dependent killing by granulocytes and monocytes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CD11a in THP-1 lysates using Integrin alpha L antibody.

注意事项 (Note)

For research use only .