

产品名称: CD127 Rabbit Monoclonal Antibody
产品货号: AMRe85393



产品概述 (Summary)

产品名称 (Production Name)	CD127 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD127
别名 (Alternative Names)	IL7R; Interleukin-7 receptor subunit alpha; IL-7 receptor subunit alpha; IL-7R subunit alpha; IL-7R-alpha; IL-7RA; CDw127; CD antigen CD127
基因 ID (Gene ID)	3575.0
蛋白 ID (SwissProt ID)	P16871.A synthetic peptide of human CD127

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 70 kDa

研究背景 (Background)

产品名称: CD127 Rabbit Monoclonal Antibody
产品货号: AMRe85393




Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway, Jak-STAT signaling pathway

图片 (Image Data)



Western blot analysis of CD127 in K562 lysates using CD127 antibody.

注意事项 (Note)

For research use only .