

产品名称: CD226 Rabbit Monoclonal Antibody
产品货号: AMRe85399



产品概述 (Summary)

产品名称 (Production Name)	CD226 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD226
别名 (Alternative Names)	Platelet and T-cell activation antigen 1; CD226
基因 ID (Gene ID)	225825.0
蛋白 ID (SwissProt ID)	Q8K4F0.Recombinant protein of mouse CD226

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 55-65 kDa

研究背景 (Background)

Involved in intercellular adhesion, lymphocyte signaling, cytotoxicity and lymphokine secretion mediated by cytotoxic T-


产品名称: CD226 Rabbit Monoclonal Antibody
产品货号: AMRe85399



lymphocyte (CTL) and NK cell. Cell surface receptor for NECTIN2. Upon ligand binding, stimulates T-cell proliferation and cytokine production, including that of IL2, IL5, IL10, IL13, and IFNG. Competes with PVRIG for NECTIN2-binding.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CD226 in mouse spleen lysates using CD226 antibody.

注意事项 (Note)

For research use only .