

产品名称: CD27 Rabbit Monoclonal Antibody
产品货号: AMRe85400



产品概述 (Summary)

产品名称 (Production Name)	CD27 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD27
别名 (Alternative Names)	CD27; CD27L receptor; LPFS2; S152; T14; TNFRSF7; TNFSF7; Tp55
基因 ID (Gene ID)	939.0
蛋白 ID (SwissProt ID)	P26842.A synthetic peptide of human CD27

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 49 kDa

研究背景 (Background)

Receptor for CD70/CD27L. May play a role in survival of activated T-cells. May play a role in apoptosis through association

产品名称: CD27 Rabbit Monoclonal Antibody
产品货号: AMRe85400




with SIVA1.

研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of CD27 in Raji, Ramos lysates using CD27 antibody.

注意事项 (Note)

For research use only .