

产品名称: CD7 Rabbit Monoclonal Antibody
产品货号: AMRe85411



产品概述 (Summary)

产品名称 (Production Name)	CD7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD7
别名 (Alternative Names)	CD7; T-cell antigen CD7; GP40; T-cell leukemia antigen; T-cell surface antigen Leu-9; TP41; CD7
基因 ID (Gene ID)	924.0
蛋白 ID (SwissProt ID)	P09564.A synthetic peptide of human CD7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 35-40 kDa

研究背景 (Background)


产品名称: CD7 Rabbit Monoclonal Antibody
产品货号: AMRe85411





Not yet known.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of CD7 in Jurkat lysates using CD7 antibody.



Immunocytochemistry analysis of CD7 (green) in Jurkat using CD7 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: CD7 Rabbit Monoclonal Antibody
产品货号: AMRe85411



注意事项 (Note)

For research use only .