

产品名称: CDA Rabbit Monoclonal Antibody
产品货号: AMRe85413



产品概述 (Summary)

产品名称 (Production Name)	CDA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDA
别名 (Alternative Names)	cytidine deaminase; CDD
基因 ID (Gene ID)	978.0
蛋白 ID (SwissProt ID)	P32320.Recombinant protein of human CDA

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa

研究背景 (Background)


This enzyme scavenges exogenous and endogenous cytidine and 2'-deoxycytidine for UMP synthesis.

产品名称: CDA Rabbit Monoclonal Antibody
产品货号: AMRe85413




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CDA in HeLa lysates using CDA antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using CDA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .