

产品名称: CDC25 Rabbit Monoclonal Antibody
产品货号: AMRe85414



产品概述 (Summary)

产品名称 (Production Name)	CDC25 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDC25
别名 (Alternative Names)	GNRP; GRF1; CDC25; GRF55; CDC25L; H-GRF55; PP13187; ras-GRF1
基因 ID (Gene ID)	5923.0
蛋白 ID (SwissProt ID)	Q13972.A synthetic peptide of human CDC25

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 145 kDa; Observed MW: 145 kDa

研究背景 (Background)

Promotes the exchange of Ras-bound GDP by GTP.


产品名称: CDC25 Rabbit Monoclonal Antibody
产品货号: AMRe85414



研究领域 (Research Area)

MAPK signaling pathway

图片 (Image Data)



Western blot analysis of CDC25 in rat Brain lysates using CDC25 antibody.

注意事项 (Note)

For research use only .