

产品名称: CECR5 Rabbit Monoclonal Antibody
产品货号: AMRe85425



产品概述 (Summary)

产品名称 (Production Name)	CECR5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CECR5
别名 (Alternative Names)	haloacid dehalogenase like hydrolase domain containing 5; CECR5
基因 ID (Gene ID)	27440.0
蛋白 ID (SwissProt ID)	Q9BXW7.A synthetic peptide of human CECR5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa; Observed MW: 46 kDa

研究背景 (Background)

MiscellaneousCandidate gene for the Cat Eye Syndrome (CES), a developmental disorder associated with the duplication of


产品名称: CECR5 Rabbit Monoclonal Antibody
产品货号: AMRe85425



a 2 Mb region of 22q11.2. Duplication usually takes in the form of a supernumerary bisatellited isodicentric chromosome, resulting in four copies of the region (represents an inv dup(22)(q11)). CES is characterized clinically by the combination of coloboma of the iris and anal atresia with fistula, downslanting palpebral fissures, preauricular tags and/or pits, frequent occurrence of heart and renal malformations, and normal or near-normal mental development.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CECR5 in K562, HeLa lysates using CECR5 antibody.

注意事项 (Note)

For research use only .