

产品名称: Cellular Apoptosis Susceptibility Rabbit Monoclonal Antibody
产品货号: AMRe85427



产品概述 (Summary)

产品名称 (Production Name)	Cellular Apoptosis Susceptibility Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cellular Apoptosis Susceptibility
别名 (Alternative Names)	CSE1L; CAS; CSE1; Exp2; Importin-alpha re-exporter; Exportin-2; XPO2
基因 ID (Gene ID)	1434.0
蛋白 ID (SwissProt ID)	P55060.A synthetic peptide of human Cellular Apoptosis Susceptibility

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 110 kDa; Observed MW: 110 kDa

研究背景 (Background)


**产品名称: Cellular Apoptosis Susceptibility Rabbit
Monoclonal Antibody**
产品货号: AMRe85427



Export receptor for importin-alpha. Mediates importin-alpha re-export from the nucleus to the cytoplasm after import substrates (cargos) have been released into the nucleoplasm. In the nucleus binds cooperatively to importin-alpha and to the GTPase Ran in its active GTP-bound form.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ular Apoptosis Susceptibility in 3T3, HeLa lysates using ular Apoptosis Susceptibility antibody.

注意事项 (Note)

For research use only .