

产品名称: CLASP1 Rabbit Monoclonal Antibody
产品货号: AMRe85437



产品概述 (Summary)

产品名称 (Production Name)	CLASP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLASP1 CLASP1; KIAA0622; MAST1; CLIP-associating protein 1; Cytoplasmic linker-associated protein 1; Multiple asters homolog 1; Protein Orbit homolog 1; hOrbit1
别名 (Alternative Names)	
基因 ID (Gene ID)	23332.0
蛋白 ID (SwissProt ID)	Q7Z460.A synthetic peptide of human CLASP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 169 kDa; Observed MW: 160 kDa

产品名称: CLASP1 Rabbit Monoclonal Antibody
产品货号: AMRe85437




研究背景 (Background)

Microtubule plus-end tracking protein that promotes the stabilization of dynamic microtubules. Involved in the nucleation of noncentrosomal microtubules originating from the trans-Golgi network (TGN). Required for the polarization of the cytoplasmic microtubule arrays in migrating cells towards the leading edge of the cell.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CLASP1 in HeLa lysates using CLASP1 antibody.

注意事项 (Note)

For research use only .