

产品名称: CLPP Rabbit Monoclonal Antibody
产品货号: AMRe85440



产品概述 (Summary)

产品名称 (Production Name)	CLPP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLPP
别名 (Alternative Names)	Endopeptidase Clp
基因 ID (Gene ID)	8192.0
蛋白 ID (SwissProt ID)	Q16740.A synthetic peptide of human CLPP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa; Observed MW: 30 kDa

研究背景 (Background)

Clp cleaves peptides in various proteins in a process that requires ATP hydrolysis. Clp may be responsible for a fairly


产品名称: CLPP Rabbit Monoclonal Antibody
产品货号: AMRe85440



general and central housekeeping function rather than for the degradation of specific substrates.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CLPP in K562, rat Brain, 3T3, HeLa lysates using CLPP antibody.

注意事项 (Note)

For research use only .