

产品名称: Coilin Rabbit Monoclonal Antibody
产品货号: AMRe85444



产品概述 (Summary)

产品名称 (Production Name)	Coilin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Coilin
别名 (Alternative Names)	CLN80; COIL; Coilin; Coilin p80
基因 ID (Gene ID)	8161.0
蛋白 ID (SwissProt ID)	P38432.A synthetic peptide of human Coilin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 63 kDa; Observed MW: 80 kDa

研究背景 (Background)

Is a component of the nuclear coiled bodies (CBS) which are involved in the function or assembly/disassembly of


产品名称: Coilin Rabbit Monoclonal Antibody
产品货号: AMRe85444



nucleoplasmic snRNPs. During mitosis, CBS disassemble, coinciding with a mitotic-specific phosphorylation of p80 coilin.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Coilin in K562 lysates using Coilin antibody.

注意事项 (Note)

For research use only .