

产品名称: Collagen VI alpha 1 Rabbit Monoclonal Antibody
产品货号: AMRe85445



产品概述 (Summary)

产品名称 (Production Name)	Collagen VI alpha 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Collagen VI alpha 1
别名 (Alternative Names)	OPLL; BTHLM1; UCHMD1
基因 ID (Gene ID)	1291.0
蛋白 ID (SwissProt ID)	P12109.Recombinant protein of human Collagen VI

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 109 kDa; Observed MW: 147 kDa

研究背景 (Background)

Collagen VI acts as a cell-binding protein.


产品名称: Collagen VI alpha 1 Rabbit Monoclonal Antibody
产品货号: AMRe85445



研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Western blot analysis of Collagen VI in Human placenta lysates using Collagen VI alpha 1 antibody.

注意事项 (Note)

For research use only .