

**产品名称:** Creatine Kinase B type Rabbit Monoclonal Antibody  
**产品货号:** AMRe85454



## 产品概述 (Summary)

|                        |                                                   |
|------------------------|---------------------------------------------------|
| 产品名称 (Production Name) | Creatine Kinase B type Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody            |
| 宿主 (Host)              | Rabbit                                            |
| 应用 (Application)       | WB,IHC,IP                                         |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                                   |

## 产品性能 (Performance)

|                     |                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                               |
| 修饰 (Modification)   | Unmodified                                                                                 |
| 同种型 (Isotype)       | IgG                                                                                        |
| 克隆 (Clonality)      | Monoclonal                                                                                 |
| 形式 (Form)           | Liquid                                                                                     |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.   |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                      |

## 免疫原信息 (Immunogen)

|                        |                                                            |
|------------------------|------------------------------------------------------------|
| 基因名 (Gene Name)        | Creatine Kinase B type                                     |
| 别名 (Alternative Names) | BCK; Ckb; CKBB; HEL 211; HEL S 29                          |
| 基因 ID (Gene ID)        | 1152.0                                                     |
| 蛋白 ID (SwissProt ID)   | P12277.A synthetic peptide of human Creatine kinase B type |

## 产品应用 (Application)

|                          |                                             |
|--------------------------|---------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 43 kDa; Observed MW: 43 kDa  |

## 研究背景 (Background)


**产品名称: Creatine Kinase B type Rabbit Monoclonal Antibody**  
**产品货号: AMRe85454**




Reversibly catalyzes the transfer of phosphate between ATP and various phosphogens (e.g. creatine phosphate). Creatine kinase isoenzymes play a central role in energy transduction in tissues with large, fluctuating energy demands, such as skeletal muscle, heart, brain and spermatozoa.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of Creatine kinase B type in K562, rat Brain, C6, 3T3 lysates using Creatine kinase B type antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon using Creatine Kinase B type antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## 注意事项 (Note)

For research use only .