

产品名称: CRM1 Rabbit Monoclonal Antibody
产品货号: AMRe85456



产品概述 (Summary)

产品名称 (Production Name)	CRM1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRM1
别名 (Alternative Names)	XPO1; CRM1; Exportin-1; Exp1; Chromosome region maintenance 1 protein homolog
基因 ID (Gene ID)	7514.0
蛋白 ID (SwissProt ID)	O14980.A synthetic peptide of human CRM1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 123 kDa; Observed MW: 123 kDa

研究背景 (Background)


产品名称: CRM1 Rabbit Monoclonal Antibody
产品货号: AMRe85456



Mediates the nuclear export of cellular proteins (cargos) bearing a leucine-rich nuclear export signal (NES) and of RNAs. In the nucleus, in association with RANBP3, binds cooperatively to the NES on its target protein and to the GTPase RAN in its active GTP-bound form (Ran-GTP). Docking of this complex to the nuclear pore complex (NPC) is mediated through binding to nucleoporins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CRM1 in 3T3 lysates using CRM1 antibody.

注意事项 (Note)

For research use only .