

产品名称: Cullin 1 Rabbit Monoclonal Antibody
产品货号: AMRe85463



产品概述 (Summary)

产品名称 (Production Name)	Cullin 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cullin 1
别名 (Alternative Names)	CUL1; Cullin1
基因 ID (Gene ID)	8454.0
蛋白 ID (SwissProt ID)	Q13616.A synthetic peptide of human Cullin 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 90 kDa; Observed MW: 90 kDa

研究背景 (Background)

Core component of multiple cullin-RING-based SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complexes,


产品名称: Cullin 1 Rabbit Monoclonal Antibody
产品货号: AMRe85463




which mediate the ubiquitination of proteins involved in cell cycle progression, signal transduction and transcription.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Cullin 1/CUL1 in MEF, 3T3 lysates using Cullin 1 antibody.



Western blot analysis of Cullin 1/CUL1 in K562, C6 lysates using Cullin 1/CUL1 antibody.



Immunocytochemistry analysis of Cullin 1(green) in HeLa using Cullin 1 antibody, and DAPI(blue)

产品名称: Cullin 1 Rabbit Monoclonal Antibody
产品货号: AMRe85463



注意事项 (Note)

For research use only .