

产品名称: Cyclin T1 Rabbit Monoclonal Antibody
产品货号: AMRe85475



产品概述 (Summary)

产品名称 (Production Name)	Cyclin T1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cyclin T1
别名 (Alternative Names)	CCNT; CYCT1; HIVE1
基因 ID (Gene ID)	904.0
蛋白 ID (SwissProt ID)	O60563.Recombinant protein of human Cyclin T1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 81 kDa; Observed MW: 81 kDa

研究背景 (Background)

Regulatory subunit of the cyclin-dependent kinase pair (CDK9/cyclin-T1) complex, also called positive transcription


产品名称: Cyclin T1 Rabbit Monoclonal Antibody
产品货号: AMRe85475




elongation factor B (P-TEFb), which is proposed to facilitate the transition from abortive to productive elongation by phosphorylating the CTD (carboxy-terminal domain) of the large subunit of RNA polymerase II (RNA Pol II).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cyclin T1 in mouse brain, rat brain lysates using Cyclin T1 antibody.



Western blot analysis of Cyclin T1 in HeLa, Jurkat, K562 lysates using Cyclin T1 antibody.

注意事项 (Note)

For research use only .