

**产品名称: Cytochrome C Oxidase subunit Vic Rabbit
Monoclonal Antibody**
产品货号: AMRe85480



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome C Oxidase subunit Vic Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytochrome C Oxidase subunit Vic
别名 (Alternative Names)	cytochrome c oxidase subunit 6C
基因 ID (Gene ID)	1345.0
蛋白 ID (SwissProt ID)	P09669.A synthetic peptide of human Cytochrome C Oxidase subunit Vic

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 9 kDa; Observed MW: 12 kDa

研究背景 (Background)


产品名称: Cytochrome C Oxidase subunit Vic Rabbit
Monoclonal Antibody
产品货号: AMRe85480




This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Cytochrome C Oxidase subunit Vic in rat heart lysates using Cytochrome C Oxidase subunit Vic antibody.



Western blot analysis of Cytochrome C Oxidase subunit Vic in 3T3, HeLa lysates using Cytochrome C Oxidase subunit Vic antibody.

**产品名称: Cytochrome C Oxidase subunit Vic Rabbit
Monoclonal Antibody**
产品货号: AMRe85480



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Cytochrome C Oxidase subunit Vic antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .