

产品名称: Cytochrome P450 Reductase Rabbit Monoclonal Antibody
产品货号: AMRe85482



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 Reductase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytochrome P450 Reductase
别名 (Alternative Names)	CPR; CYPOR; P450R; por
基因 ID (Gene ID)	5447.0
蛋白 ID (SwissProt ID)	P16435.A synthetic peptide of human Cytochrome P450 Reductase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 77 kDa


研究背景 (Background)

产品名称: Cytochrome P450 Reductase Rabbit Monoclonal Antibody
产品货号: AMRe85482


This enzyme is required for electron transfer from NADP to cytochrome P450 in microsomes. It can also provide electron transfer to heme oxygenase and cytochrome B5.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Cytochrome P450 Reductase in K562, rat Brain, C6, 3T3 lysates using Cytochrome P450 Reductase antibody.



Immunocytochemistry analysis of Cytochrome P450 Reductase (green) in HeLa using Cytochrome P450 Reductase antibody, and DAPI (blue)

产品名称: Cytochrome P450 Reductase Rabbit Monoclonal
Antibody
产品货号: AMRe85482



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Cytochrome P450 Reductase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .