

产品名称: DCAMKL1 Rabbit Monoclonal Antibody
产品货号: AMRe85494



产品概述 (Summary)

产品名称 (Production Name)	DCAMKL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DCAMKL1
别名 (Alternative Names)	doublecortin-like kinase 1; CL1; DCLK; CLICK1; DCDC3A; DCAMKL1
基因 ID (Gene ID)	9201.0
蛋白 ID (SwissProt ID)	O15075.A synthetic peptide of human DCAMKL1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 82 kDa; Observed MW: 82 kDa

研究背景 (Background)

Probable kinase that may be involved in a calcium-signaling pathway controlling neuronal migration in the developing


产品名称: DCAMKL1 Rabbit Monoclonal Antibody
产品货号: AMRe85494




brain. May also participate in functions of the mature nervous system.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of DCAMKL1 in mouse brain lysates using DCAMKL1 antibody.




Western blot analysis of DCAMKL1 in rat Brain lysates using DCAMKL1 antibody



Immunocytochemistry analysis of DCAMKL1 (green) in U87-MG using DCAMKL1 antibody, and DAPI(blue).

产品名称: DCAMKL1 Rabbit Monoclonal Antibody
产品货号: AMRe85494



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using DCAMKL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .