

产品名称: DcR2 Rabbit Monoclonal Antibody
产品货号: AMRe85496



产品概述 (Summary)

产品名称 (Production Name)	DcR2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DcR2 TNFRSF10D; DCR2; TRAILR4; TRUNDD; Tumor necrosis factor receptor superfamily member 10D; Decoy receptor 2; DcR2; TNF-related apoptosis-inducing ligand receptor 4; TRAIL receptor 4; TRAIL-R4; TRAIL receptor with a truncated death domain; CD
别名 (Alternative Names)	
基因 ID (Gene ID)	8793.0
蛋白 ID (SwissProt ID)	Q9UBN6.A synthetic peptide of human DcR2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 35 kDa

产品名称: DcR2 Rabbit Monoclonal Antibody
产品货号: AMRe85496




研究背景 (Background)

Receptor for the cytotoxic ligand TRAIL. Contains a truncated death domain and hence is not capable of inducing apoptosis but protects against TRAIL-mediated apoptosis. Reports are contradictory with regards to its ability to induce the NF-kappa-B pathway.


研究领域 (Research Area)

Apoptosis, TGF-beta signaling pathway

图片 (Image Data)



Western blot analysis of DcR2 in Jurkat lysates using DcR2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DcR2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .